


# Refurbished Lyotherm2 Available



This system has been fully upgraded to current manufacturing standards, and is supplied fully tested and validated ensuring the system is operating as per current specifications

## Included in this offer:

- Refurbished Lyotherm2 Control Unit
- Refurbished Thermal Heating Block
- New PT100 Temperature Probe Set
- New Impedance Probe
- Microsoft Office 2016
- New PC (i5, 8GB RAM, Windows 10)



DTA: Onset of crystallisation at -54.1°C, second crystallisation at -26°C, Z: Tg' at -31.3°C

Only one system is for sale and is offered subject to availability. Shipment and installation is included.

**For more information, please contact:**  
**tcodd@biopharma.co.uk | +44 (0) 1962 841092**